



Government of India

NATIONAL CENTRE FOR DISEASE CONTROL

(DIRECTORATE GENERAL OF HEALTH SERVICES)

22, SHAM NATH MARG, NEW DELHI – 110054

CENTRE FOR ARBOVIRAL AND ZOONOTIC DISEASES



Proforma for sample submission of a suspected case of Echinococcus/ Hydatid Disease for serology testing

Lab Reference No:

Date:

Hospital:

Hospital IP/OP number:

Name:

Father's/Husband's Name:

Age/Gender:

Occupation:

Complete address:

Telephone/ Mobile number:

E-mail ID:

HISTORY OF PRESENT ILLNESS

Asymptomatic: ☐ Yes/ ☐ No

Fever: ☐ Present/ ☐ Absent

Date of onset of fever:

Course of fever: ☐ Intermittent/ ☐ Remittent/ ☐ Continuous

Type of fever: ☐ Low grade/ ☐ High grade

Abdominal Pain: ☐ Present/ ☐ Absent

Site:

Duration:

Tenderness: ☐ Present/ ☐ Absent

Duration:

Abdominal Swelling: ☐ Present/ ☐ Absent

Site:

Duration:

Jaundice: ☐ Present/ ☐ Absent

Duration:

Nausea/Vomiting: ☐ Present/ ☐ Absent

Hepatomegaly ☐ Present/ ☐ Absent

Splenomegaly

☐ Present/ ☐ Absent

Pressure symptoms due to swelling: ☐ Present/ ☐ Absent

Any other symptoms (extra hepatic):

HISTORY OF PAST ILLNESS

History of contact with dog: ☐ Present/ ☐ Absent

Duration:

FAMILY HISTORY

Any other person in the family suffering from the similar illness?

☐ Yes/ ☐ No

Detailed medical history (if yes)

LABORATORY INVESTIGATION

Eosinophil count:

Serum bilirubin:

SGOT:

SGPT:

Findings of Plain X-Ray Abdomen (attach report):

Findings of Ultrasound of Abdomen (attach report):

Findings of CT scan if done (attach report)

Clinical Finding:

Date of sample collection*:

Signature, Name, Designation, and contact number of consultant doctor with Seal

***Guidelines for sample collection and transportation**

SAMPLE: 1-2 mL non-hemolyzed serum in sterile plain vial/ 5 mL blood in clot activator and gel vial/ 5 mL blood in vial kept for 30-40 mins for clot retraction.

PACKAGING AND TRANSPORTATION: For serum, triple layer package to be sent in cold chain maintained at 2-8°C for transportation times less than 48 hours and at -20°C for transportation times more than 48 hours.